

BP5

- 5% Phenyl methylpolysiloxane

ID(mm)	Length	Film	Iso/Prog	Cat. No.	
0.22	12m	0.25µm	340/350	054167	
	25m	0.25µm	340/350	054168	
	30m	0.25µm	340/350	054171	
	50m	0.25µm	340/350	054169	
	50m	1.00µm	340/350	054175	
0.25	15m	0.25µm	340/350	054182	
	30m	0.25µm	340/350	054183	
	30m	1.00µm	340/350	054203	
	60m	0.25µm	340/350	054184	
	60m	1.00µm	340/350	054215	
0.32	12m	0.25µm	340/350	054179	
	15m	0.25µm	340/350	054176	
	25m	0.25µm	340/350	054180	
	25m	0.50µm	340/350	054186	
	25m	1.00µm	340/350	054192	
	30m	0.25µm	340/350	054177	
	30m	0.50µm	340/350	054216	
	30m	1.00µm	340/350	054189	
	50m	0.50µm	340/350	054187	
	50m	1.00µm	340/350	054193	
	60m	0.25µm	340/350	054178	
	60m	1.00µm	340/350	054188	
	0.53	12m	1.00µm	340/350	054197
		15m	1.00µm	340/350	054194
		15m	1.50µm	280/300	054199
25m		1.00µm	340/350	054198	
30m		0.50µm	340/350	0541935	
30m		1.00µm	340/350	054195	
30m		5.00µm	280/300	054196	
60m	1.50µm	340/350	054204		

BPX5

- 5% Phenyl polysilphenylene-siloxane
- Ultra low bleed

ID(mm)	Length	Film	Iso/Prog	Cat. No.	
0.10	10m	0.10µm	360/370	054099	
	0.15	10m	1.20µm	054106	
0.15	12m	0.25µm	360/370	054103	
	12m	0.40µm	360/370	054107	
	25m	0.25µm	360/370	054104	
	25m	0.40µm	360/370	054108	
	30m	0.15µm	360/370	054110	
	50m	0.25µm	360/370	054105	
	0.22	12m	0.25µm	360/370	054112
		25m	0.25µm	360/370	054113
		25m	1.00µm	360/370	054116
		30m	0.25µm	360/370	054142
50m		0.25µm	360/370	054114	
50m		1.00µm	360/370	054117	
0.25		15m	0.25µm	360/370	054100
	15m	1.00µm	360/370	054121	
	30m	0.25µm	360/370	054101	
	30m	0.50µm	360/370	0541025	
	30m	1.00µm	360/370	054122	
	60m	0.25µm	360/370	054102	
	60m	1.00µm	360/370	054123	
	0.32	6m	1.00µm	360/370	0541261
		12m	0.25µm	360/370	054118
		12m	0.50µm	360/370	054124
12m		1.00µm	360/370	054127	
15m		0.25µm	360/370	054144	
15m		1.00µm	360/370	054152	
25m		0.25µm	360/370	054119	
25m		0.35µm	360/370	054125	
25m		1.00µm	360/370	054128	
25m		3.00µm	350/360	054136	
30m		0.25µm	360/370	054145	
30m		0.50µm	360/370	0541205	
30m		1.00µm	360/370	054153	
50m		0.25µm	360/370	054120	
50m		0.50µm	360/370	054126	
50m	1.00µm	360/370	054129		
60m	0.25µm	360/370	054146		
60m	1.00µm	360/370	054154		
0.53	12m	0.25µm	360/370	054133	
	12m	1.00µm	360/370	054130	
	12m	3.00µm	350/360	054138	
	15m	0.50µm	360/370	0541344	
	15m	1.00µm	360/370	054147	
	15m	1.50µm	350/360	0541347	
	15m	3.00µm	350/360	054159	
	25m	0.25µm	360/370	054134	
	25m	1.00µm	360/370	054131	
	25m	3.00µm	350/360	054139	
	30m	0.50µm	360/370	0541345	
	30m	1.00µm	360/370	054148	
	30m	1.50µm	350/360	0541348	
	30m	3.00µm	350/360	054160	
	50m	1.00µm	360/370	054132	
60m	1.00µm	360/370	054158		

