

















GC Supplies




SGE Deactivated Inlet Liners

Description	(mm) ID	(mm) Length	5/pk Cat. No.	25/pk Cat. No.
HP5890, HP6890, HP6850 & HP4890				
 Split / Splitless FocusLiner	4.0	78.5	092002	092219
 Split / Splitless Tapered FocusLiner	4.0	78.5	092003	092011
 Split / Splitless Fast FocusLiner	2.3	78.5	092005	092008
 Split / Splitless Tapered Fast FocusLiner	2.3	78.5	092111	092115
 Split, straight-through liner	4	78.5	092007	092222
 Split with quartz wool	4	78.5	092001	092220
 Split / Splitless with single taper	4	78.5	092017	092229
 Split / Splitless with single taper and quartz wool	4	78.5	092019	
 Split / Splitless with double taper	4	78.5	092018	092230
 Direct, straight-through liner	1.2	78.5	092016	092224
 Splitless (quartz) straight-through liner	2	78.5	092004	
 Splitless with recessed gooseneck	2	78.5	092013	
 Recessed gooseneck with quartz wool	4	78.5	092010	092223
 Packed Column Liner for HP5880, 5890, 6890	1.8	92	092234	092235





Perkin Elmer AutoSystem

 Split / Splitless FocusLiner	4	92	092092	
 Split / Splitless Tapered FocusLiner	4	92	092095	






Shimadzu 17A (SPL-17 injector)

 Split / Splitless FocusLiner	3.4	95	092062	
 Split / Splitless Tapered FocusLiner	3.4	95	092068	
 Split / Splitless with middle gooseneck	3.4	95	092085	




Shimadzu 14/15A/16 (SPL-14 injector)

 Split / Splitless FocusLiner	3.4	99	092065	
 Split / Splitless Tapered FocusLiner	3.4	99	092066	
 Split / Splitless with 2.0mm middle gooseneck	3.4	99	092082	
 Split / Splitless with single taper	3.4	99	0920831	



ThermoQuest - CE Instruments Model 8000/TRACE

 Split / Splitless FocusLiner	5	105	092045	
 Split / Splitless, Tapered FocusLiner	5	105	092046	
 Splitless with single taper	3	105	092141	
 Split with single taper	5	105	092144	
 Splitless straight-through liner	3	105	092147	

Varian Injector Models 1075/1077

 Split / Splitless FocusLiner	4	72	092022	
 Split / Splitless Tapered FocusLiner	4	72	092025	
 Split with quartz wool	4	72	092021	

Varian Injector Models 1078/1079

 Split / Splitless FocusLiner	3.4	54	092037	
 Split / Splitless Tapered FocusLiner	3.4	54	092036	