

GC Supplies



Column Install Kit



Column Rinse Kit

GC Column Installation Kit

The capillary column installation kit contains:

- Diamond-tipped pencil
- 20x magnifier
- Ferrule removal tool
- Septum removal tool
- Refillable 2 oz bottle of Snoop

Cat. No.	Description
430-2000	Kit with diamond-tipped cutting pencil

GC Column Rinse Kit

This kit makes solvent rinsing easy! Solvent rinsing removes many contaminants that accumulate in columns. The contaminants are removed by passing solvents through the column. Only bonded and cross-linked phase columns should be rinsed, or permanent damage to the column will result.

Cat. No.	Description
430-3000	Capillary rinse kit

Injection Port Cleaning Kit



- 4 Stainless Steel brushes: 1/8", 3/16", 1/4" and 3/8"
- Scraper to remove residue

Cat. No.	Description
124500	Injector cleaning kit

FID Cleaning Kit



This easy-to-use maintenance kit improves sensitivity, reduces noise and decreases detector spiking.

Cat. No.	Description
9301-0985	FID cleaning kit

LD-228 Helium Leak Detector

- Detection limits as low as 0.01mL/min. for He



This rechargeable gas leak detector is designed for the detection of helium gas. It is fast, reliable and very easy to use. The pen-sized sample probe proves handy in confined spaces.

The detection principle of the new LD-228 Helium Gas Leak Detector is based on the difference in thermal conductivity between the sample gas and atmosphere. The detector contains a high sensitivity thermistor actuated cell and a vacuum sample intake.

When a leak is detected through the monitoring probe, an elevated reading occurs on the LED level meter and an alarm buzzer sounds.

Cat. No.	Description
AL60229	LD-228 Helium gas leak detector

Snoop

Leak test most pressurized gas systems. We recommend that leak detecting solutions be used with care where there is any possibility of the fluid being aspirated into gas flow streams supplying the column or detector. This squeezable bottle comes with a "Snooper Tube" extendable up to 12 inches.

Cat. No.	Description
9300-0311	Snoop - 8 ounce bottle

Glas-Treet

- Used to deactivate active sites found in glassware, chromatography columns, and glass wool
- Does not coat glass; renders it inactive chemically

Cat. No.	Description
970069	Glas-Treet, 4 x 25mL

Col-Treet

- Silylating mixture designed for conditioning gas chromatography columns
- Eliminates peak tailing, peak broadening and residues which may reduce column efficiency

Col-Treet is NOT to be used with any hydrogen bonding stationary phases such as CARBOWAX® or phases with active hydrogen sites.

Cat. No.	Description
970130	Col-Treet, 10 x 1mL
970131	Col-Treet, 4 x 5mL