

# GC Supplies

## Agilent Instrument Parts and Supplies

### MSD Ion Source Cleaning Supplies

Cat. No.	Description	Qty.
8650-0030	Nylon gloves, lint-free, large	1
8650-0029	Nylon gloves, lint-free, small	1
9310-4828	Lint-free wipes, cotton, 9 x 9"	300
05980-60051	Cloths, lint-free	15
5061-5896	Abrasive sheets, 600 mesh	5
8500-0656	PFTBA sample, certified	10g
05980-20018	PFTBA glass vial	1
5080-5400	Cotton swabs	100

### 5972/5971 (EI) MSD Ion Source Parts


Cat. No.	Description
05971-20126	Entrance lens
05971-20130	Lens insulator, 2/pk
05971-20143	Ion focus lens
G1072-20008	Drawout cylinder
05971-20134	Drawout plate
05971-20128	Ion source body
0515-1446	Set screw
05971-60170	Repeller assembly
0515-1046	Screw (for filament on the source)
05971-20305	Transfer line tip, gold plated

### MSD Filaments

Cat. No.	Description
05972-60053	Filament, 5973 (EI), 5972 (EI/CI)
G1099-80053	Filament, 5973 (CI), 2/pk
05971-60140	Filament, 5971 (EI/CI), GCD (EI)

### 5973 (EI) MSD Ion Source Parts

Cat. No.	Description
G1099-20130	Ion source body
G1099-20132	Repeller
G1099-20136	Interface socket
05971-20134	Drawout plate
G1072-20008	Drawout cylinder
05971-20143	Ion focus lens
05971-20126	Entrance lens
G1099-20133	Repeller insulator
0515-1446	Set screw
0535-0071	Repeller nut, 5.5mm, ultraclean
G1099-60104	Ion source sensor
05971-20130	Lens insulator, 2/pk
G1099-60106	Ion source assembly
G1099-60177	Source heater assembly
G1099-60170	Repeller assembly



5973 Ion Source Assembly (EI)



Cleaning and Maintenance Supplies



5973/5972 Filament

### MSD Vacuum System and Pumps

#### Calibration samples

Cat. No.	Description
8500-0656	PFTBA sample, certified, 10 grams
8500-8130	PFDTD sample, 10 g, for 5973 CI source

#### Diffusion pump

6040-0809	Fluid, SantoVac Ultra 5P, 18.5mL
0960-0376	Ion gauge tube for measuring vacuum (71/72)
0960-0897	Triode gauge tube for meas. vacuum (72/73)

#### Foreline pump

6040-0798	Rough pump oil, 1 gal Inland 45
6040-0834	Rough pump oil, 1 quart Inland 45
9301-1104	Molecular sieve



### Electron Multiplier Replacement Horn

Cat. No.	Description
05971-80103	Electron multiplier replacement horn