

GC Supplies

SGE Instrument Supplies

Agilent 5890 and 6890



Autosampler Syringes, Fixed Needle for 7673 and 7683

Cat. No.	Volume	Length	Gauge	OD	Tip
002810	10µL	42	23	0.63	Cone
002821	10µL	42	23-26	0.63/0.47	Cone
002822	10µL, 6/pk	42	23-26	0.63/0.47	Cone

Septa

Cat. No.	Description	Qty.
041882	Auto-Sep T, 11mm	25



Liners for Split/Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092002	FocusLiner	4.0	6.3	78.5	5
092003	FocusLiner, tapered	4.0	6.3	78.5	5
092017	Single Taper Liner	4.0	6.3	78.5	5
092007	Split, Straight Liner	4.0	6.3	78.5	5
092222	Split, Straight Liner	4.0	6.3	78.5	25
0726532	Viton O-Ring	-	-	-	10

15% Graphite/85% Vespel Ferrules (not for GC/MS)

Cat. No.	Ferrule ID	Column ID	Qty.
073109	0.4mm	.10 - .25	10
073111	0.5mm	.32	10
073113	0.8mm	.53	10

15% Graphite/85% Vespel Ferrules (for GC/MS)

072663	0.4mm	.1-.25	10
--------	-------	--------	----

100% Graphite Ferrules (not for GC/MS)

072635	0.5mm	.1-.32	10
--------	-------	--------	----




FID Flame Jets

Cat. No.	Description
103474	FID flame jet for Agilent 5890, 61.3mm
103475	FID flame jet for Agilent 6890, 42.3mm

Electron Multipliers

Cat. No.	MS Description
14511	5970
14516	5971/5972
14616	5973 (multiplier only)
14617	5973 (first time installation – includes mount)



Perkin Elmer AutoSystem and Clarus 500

Autosampler Syringes, Fixed Needle

Cat. No.	Volume	Length	Gauge	OD	Tip
001957	5µL	70	23	0.63	Cone

Septa

Cat. No.	Description	Qty.
041882	Auto-Sep T, 11mm	25

Liners for Split/Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092092	FocusLiner	4.0	6.2	92	5
092095	FocusLiner, tapered	4.0	6.2	92	5

15% Graphite/85% Vespel Ferrules (not for GC/MS)

Cat. No.	Ferrule ID	Column ID	Qty.
072663	0.4mm	.10 - .25	10
072654	0.5mm	.32	10
072655	0.8mm	.53	10

Electron Multipliers

Cat. No.	MS Description
14532	Q-Mass