

















PerkinElmer Brownlee™ Cartridge Holders

Ordering Your First Cartridge Column?

If you are ordering a cartridge for the first time, you also need to order a holder of the same length. A holder includes all parts necessary for use: two end fittings, body, and two sets of nuts and ferrules.

The Universal Holder Kit is a cost effective way of having one complete holder or system holder of any size with a coupled guard cartridge.



	Cat. No.	Description
	0715-0025	Universal Holder Kit, Includes: 3, 10, and 22cm Bodies; 2 End assemblies, 1 union, and 1 NewGuard end assembly
	0715-0001	NewGuard™ Holder, Holds a single NewGuard cartridge for use with conventional HPLC columns
	0715-0013	3cm Holder Holds a single 3cm cartridge
	0715-0014	10cm Holder Holds a single 10cm cartridge
	0715-0015	22cm Holder Holds a single 22cm cartridge
	0715-0016	10cm NewGuard System Holder For direct coupling of NewGuard and 10cm cartridge
	0715-0017	22cm NewGuard System Holder For direct coupling of NewGuard and 22cm cartridge
	0715-0018	Union, For direct coupling of two cartridges. (2 Holder bodies and 2 end assemblies also required)
	0715-0002	NewGuard End Assembly Upgrades 3cm guard holder to NewGuard Holder
	0715-0019	End Assembly Holder component
	0715-0020	3cm Holder Body Holder component
	0715-0021	10cm Holder Body Holder component
	0715-0022	22cm Holder Body Holder component
	0715-0023	Seal Replacement Kit, Tools and two 7000 psi seals. <i>Not needed first year.</i>
	0715-0024	Two 7000 psi seals (Seal replacement kit required.) <i>Not needed first year.</i>
	0715-0004	3cm Prep-10 Holder Includes end assemblies and body
	0715-0005	10cm Prep-10 Holder Includes end assemblies and body
	0715-0006	25cm Prep-10 Holder Includes end assemblies and body
	0715-0011	Prep-10 Seal Replacement Kit, tools only, seals not included
	0715-0012	Prep-10 Inlet Seal, 1 each
	0715-0026	Prep-10 Outlet Seal, 1 each