

# GC Supplies

## SGE Instrument Supplies

### Agilent 5890 and 6890



#### Autosampler Syringes, Fixed Needle for 7673 and 7683

Cat. No.	Volume	Length	Gauge	OD	Tip
002810	10µL	42	23	0.63	Cone
002821	10µL	42	23-26	0.63/0.47	Cone
002822	10µL, 6/pk	42	23-26	0.63/0.47	Cone

#### Septa

Cat. No.	Description	Qty.
041882	Auto-Sep T, 11mm	25



#### Liners for Split/Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092002	FocusLiner	4.0	6.3	78.5	5
092003	FocusLiner, tapered	4.0	6.3	78.5	5
092017	Single Taper Liner	4.0	6.3	78.5	5
092007	Split, Straight Liner	4.0	6.3	78.5	5
092222	Split, Straight Liner	4.0	6.3	78.5	25
0726532	Viton O-Ring	-	-	-	10

#### 15% Graphite/85% Vespel Ferrules (not for GC/MS)

Cat. No.	Ferrule ID	Column ID	Qty.
073109	0.4mm	.10 - .25	10
073111	0.5mm	.32	10
073113	0.8mm	.53	10

#### 15% Graphite/85% Vespel Ferrules (for GC/MS)

072663	0.4mm	.1-.25	10
--------	-------	--------	----

#### 100% Graphite Ferrules (not for GC/MS)

072635	0.5mm	.1-.32	10
--------	-------	--------	----




#### FID Flame Jets

Cat. No.	Description
103474	FID flame jet for Agilent 5890, 61.3mm
103475	FID flame jet for Agilent 6890, 42.3mm

#### Electron Multipliers

Cat. No.	MS Description
14511	5970
14516	5971/5972
14616	5973 (multiplier only)
14617	5973 (first time installation – includes mount)



### Perkin Elmer AutoSystem and Clarus 500

#### Autosampler Syringes, Fixed Needle

Cat. No.	Volume	Length	Gauge	OD	Tip
001957	5µL	70	23	0.63	Cone

#### Septa

Cat. No.	Description	Qty.
041882	Auto-Sep T, 11mm	25

#### Liners for Split/Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092092	FocusLiner	4.0	6.2	92	5
092095	FocusLiner, tapered	4.0	6.2	92	5

#### 15% Graphite/85% Vespel Ferrules (not for GC/MS)

Cat. No.	Ferrule ID	Column ID	Qty.
072663	0.4mm	.10 - .25	10
072654	0.5mm	.32	10
072655	0.8mm	.53	10

#### Electron Multipliers

Cat. No.	MS Description
14532	Q-Mass

### Shimadzu

#### Autosampler Syringes, Removeable Needle for AOC 14, 17, 20

Cat. No.	Volume	Length	Gauge	OD	Tip
002898	10µL	42	23	0.63	Cone

#### Septa for 17A, 2010

Cat. No.	Description	Qty.
041890	Shimadzu Plug, Enduro Blue	50



#### Liners for Split/Splitless Injection for 17A, 2010

Cat. No.	Description	ID	OD	Length	Qty.
092062	FocusLiner	3.4	5.0	95	5
092068	FocusLiner, tapered	3.4	5.0	95	5

#### 100% Graphite Ferrules (not for GC/MS)

Cat. No.	Ferrule ID	Column ID	Qty.
0726080	0.5mm	.10 - .25	10
0726082	0.8mm	.45 - .53	10



#### Electron Multipliers

Cat. No.	MS Description
14533	QP5000

### Thermo Finnigan TRACE GC



#### Autosampler Syringes, Fixed Needle for AS800 and AS2000

Cat. No.	Volume	Length	Gauge	OD	Tip
001981	5µL	50	23	0.63	Cone
002980	10µL	50	26	0.47	Cone
002992	10µL	80	26	0.47	Cone

#### Septa

Cat. No.	Description	Qty.
041886	Auto-Sep T, 17mm	25



#### Liners for Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092141	Single Taper Liner	3.0	8.0	105	5

### Thermo Finnigan TRACE GC (continued)

#### SilTite Metal Ferrules

Cat. No.	Ferrule ID	Column ID	Qty.
073450	0.4mm	.10 - .25	10
073451	0.5mm	.32	10
073452	0.8	.53	10

#### Electron Multipliers

Cat. No.	Description
14139	4000 (PPINICI)
14138	Bench-top ion pre GC-Q

### Varian



#### Autosampler Syringes, Removeable Needle for 8200

Cat. No.	Volume	Length	Gauge	OD	Tip
002924	10µL	53	25	0.50	S/Hole

#### Septa for 1075, 1077, 1078, 1079, 1093, 1094

Cat. No.	Description	Qty.
041882	Auto-Sep T, 11mm	25

#### Septa for 1177

Cat. No.	Description	Qty.
041879	Auto-Sep T, 9mm	25

#### Liners for Model 1078/1079 Injector

Cat. No.	Description	ID	OD	Length	Qty.
092037	FocusLiner™	3.4	8.0	105	5

#### Liners for Split/Splitless Injection

Cat. No.	Description	ID	OD	Length	Qty.
092017	Single Taper Liner	4.0	6.3	78.5	5
092019	Single Taper, wool	4.0	6.3	78.5	5

#### 15% Graphite/85% Vespel Ferrules – for Capillary Columns

Cat. No.	Ferrule ID	Column ID	Qty.
072663	0.4mm	.10 - .25	10
072654	0.5mm	.32	10
072655	0.8mm	.53	10



#### FID Flame Jets

Cat. No.	Description
103469	FID flame jet for Varian

#### Electron Multipliers

Cat. No.	Description
14138	Saturn (Pre-2000)
14147	Saturn 200, 2100, 2200