





Transgenomic HPLC Carbohydrate Analysis Columns



CARBOsep CHO-682 Lead

- Lead form ligand-exchange
- Ideal for the separation of mono and disaccharides as well as alcohols
- High capacity

Cat. No.	Description
CHO-99-8854	CHO-682, 7.8 x 200mm
CHO-99-9854	CHO-682, 7.8 x 300mm
CHO-99-2354	CHO-682 guard kit (1 holder, 1 cartridge)
CHO-99-1354	CHO-682 guard cartridge, 2/pk

CARBOsep CHO-820 Calcium

- Calcium form ligand-exchange
- Designed with balance of resolution and ruggedness





Cat. No.	Description
CHO-99-8855	CHO-820, 7.8 x 200mm
CHO-99-9855	CHO-820, 7.8 x 300mm
CHO-99-2355	CHO-820 guard kit (1 holder, 1 cartridge)
CHO-99-1355	CHO-820 guard cartridge, 2/pk

CARBOsep CHO-620

- Calcium form ligand exchange column
- Ideal for the separation of monosaccharides and sugar alcohols
- Very reproducible

Cat. No.	Description
CHO-99-9753	CHO-620, 6.5 x 300mm
CHO-99-2353	CHO-620 guard kit (1 holder, 1 cartridge)
CHO-99-1353	CHO-620 guard cartridge, 2/pk

Transgenomic HPLC Columns



Carbohydrate Analysis Columns

CARBOSep Coregel-87C

- Calcium form 9µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87C
- Designed for the analysis of sugars and sugar alcohols

Cat. No.	Description
CHO-99-9860	Coregel-87C, 7.8 x 300mm
CHO-99-2360	87C guard kit (1 holder, 1 cartridge)
CHO-99-1360	87C guard cartridge, 2/pk

CARBOSep Coregel-87K

- Potassium form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87K
- Target application corn syrup and molasses

Cat. No.	Description
CHO-99-9862	Coregel-87K, 7.8 x 300mm
CHO-99-1362	Coregel-87K guard cartridge, 2/pk
AXC-99-1300	Universal guard cartridge holder

CARBOSep Coregel-87N

- Sodium form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87N
- Designed for the fast separation of monosaccharides and sugar alcohols

Cat. No.	Description
CHO-99-9863	Coregel-87N, 7.8 x 300mm
CHO-99-1363	Coregel-87N guard cartridge, 2/pk
AXC-99-1300	Universal guard cartridge holder

CARBOSep Coregel-87P

- Lead form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87P
- Optimized for the analysis of cellulose hydrolysates

Cat. No.	Description
CHO-99-9864	Coregel-87P, 7.8 x 300mm
CHO-99-1364	Coregel-87P guard cartridge, 2/pk
AXC-99-1300	Universal guard cartridge holder

Organic Acid Analysis Columns

- Separates organic acids, carbohydrates, and alcohols on the same column

ICSep Coregel-87H3

- Very rugged design which provides long life

Cat. No.	Description
ICE-99-9861	Coregel-87H3 6.5 x 300mm
ICE-99-2361	87H guard kit (1 holder, 1 cartridge)
ICE-99-2371	87H guard cartridge, 2/pk

ICSep ION-300

- Select when high resolution is the primary concern

Cat. No.	Description
ICE-99-9850	ION-300, 7.8 x 300mm
ICE-99-2354	GC-801 guard kit (1 holder, 1 cartridge)
ICE-99-2364	GC-801 guard cartridge, 2/pk

ICSep ORH-801

- Provides good balance of high efficiency and ruggedness

Cat. No.	Description
ICE-99-9754	ORH-801, 6.5 x 300mm
ICE-99-2354	GC-801 guard kit (1 holder, 1 cartridge)
ICE-99-2364	GC-801 guard cartridge, 2/pk