

## Barbed Adapters

### Thread to Barbed Adapters



### Swivel Barb Adapters

➔ Barb connection spins freely from the nut to prevent twist during installation



### Luer to Barbed Adapters

Barb to Slip-Type Male Luer Adapters

Barb to Male Luer with Lock Ring Adapters

Barb to Female Luer-Lock Connectors



BARBED ADAPTERS				
PART NO	TUBING ID	THREADS	THRU-HOLE	MATERIAL
<b>THREAD TO BARBED ADAPTERS</b>				
P-639	1/16" (1.55 mm)	10-32 Flat Btm	0.02" (0.5 mm)	PEEK
P-646	1/16" (1.55 mm)	1/4-28 Flat Btm	0.04" (1.0 mm)	ETFE
P-647	1/8" (3.20 mm)	1/4-28 Flat Btm	0.08" (2.0 mm)	ETFE
P-648	3/16" (4.75 mm)	1/4-28 Flat Btm	0.10" (2.5 mm)	ETFE
P-663	1/16" (1.55 mm)	10-32 Coned	0.04" (1.0 mm)	ETFE
P-664	1/8" (3.20 mm)	10-32 Coned	0.04" (1.0 mm)	ETFE
P-665	1/8" (3.20 mm)	10-32 Coned	0.04" (1.0 mm)	ETFE
P-668	1/16" (1.55 mm)	1/4-28 Flat Btm	0.04" (1.0 mm)	PEEK
P-689	1/16" (1.55 mm)	5/16-24 Flat Btm	0.10" (2.5 mm)	ETFE
P-690	1/16" – 3/32" (1.55 – 2.38 mm)	1/4-28 Flat Btm	0.04" (1.0 mm)	ETFE
P-691	3/32" – 1/8" (2.38 – 3.20 mm)	1/4-28 Flat Btm	0.06" (1.50 mm)	ETFE
P-692	0.020" – 1/32" (0.50 – 0.80 mm)	1/4-28 Flat Btm	0.02" (0.5 mm)	PEEK
<b>SWIVEL BARB ADAPTERS</b>				
D-646	1/16" (1.55 mm)	1/4-28 Flat Btm	0.03" (0.75 mm)	Polypropylene
D-647	3/32" (2.40 mm)	1/4-28 Flat Btm	0.056" (1.5 mm)	Polypropylene
D-648	1/8" (3.20 mm)	1/4-28 Flat Btm	0.08" (2.0 mm)	Polypropylene
<b>LUER TO BARB ADAPTERS</b>				
P-854	1/16" (1.55 mm)	Male Luer (Slip-type)	0.046"	Polypropylene
P-850	1/16" (1.55 mm)	Male Luer with Lock Ring	0.049"	Polypropylene
P-851	3/32" (2.40 mm)	Male Luer with Lock Ring	0.071"	Polypropylene
P-852	1/8" (3.20 mm)	Male Luer with Lock Ring	0.099"	Polypropylene
P-853	5/32" (4.0 mm)	Male Luer with Lock Ring	0.108"	Polypropylene
P-857	1/16" (1.55 mm)	Female Luer Connection	0.030"	Polypropylene
P-858	3/32" (2.40 mm)	Female Luer Connection	0.056"	Polypropylene
P-859	1/8" (3.20 mm)	Female Luer Connection	0.080"	Polypropylene
P-870	1/16" (1.55 mm)	Female Luer Connection	0.030"	ETFE
P-871	3/32" (140 mm)	Female Luer Connection	0.056"	ETFE
P-872	1/8" (3.20 mm)	Female Luer Connection	0.080"	ETFE