

Agilent Vacuum Degassers

A vacuum degasser is recommended for:

- ➔ Maximum sensitivity in the low UV wavelength range
- ➔ High injection precision
- ➔ High retention time reproducibility
- ➔ Flow rates below 0.5 mL/min

AGILENT VACUUM DEGASSERS		
PART NO	DESCRIPTION	COMMENTS
G1322-67300	Tubing kit, degasser to pump, 4/pk, 30 cm pieces of tubing with screws and bushings	For G1322A, G1379A/B
5067-5378	Inlet tubing	For Agilent 1260 Infinity Quaternary Pump VL
0100-1710	Mounting tool for flangeless nut	For G1322A, G1379A/B
8710-1930	Plastic tubing cutter	For G1322A, G1379A/B
5063-6598	Tefzel ferrules and SS lock rings, 1/8", 10/pk	For G1322A, G1379A/B
5063-6599	PPS nuts, 1/8", 1/4-28 thread, 10/pk	For G1322A, G1379A/B
5022-2155	Union, 1/4-28 thread, polypropylene	For G1322A, G1379A/B
5062-2483	PTFE solvent tubing, 5 m, 1.5 mm ID, 3 mm OD	For G1322A, G1379A/B
5190-5103	Disposable syringes, 20 mL, 100/pk	For G1322A
9301-1337	Syringe adapter, 1/16" OD, 2" long	For G1322A
G1379-67310	Micro vacuum degasser tubing kit	For G1379A
5042-8922	Micro vacuum degasser tubing kit	For G1379B
G1322-68705	Online degasser accessory kit, Includes 8 screws, 8 bushings, 4 markers, tubing, syringe and syringe adapter	For G1322A, G1379A/B
5067-4798	Vacuum chamber (2 channel)	For G1379A, G1379B
5067-5380	Tubing kit for 5067-4798 vacuum chamber	For G1379A, G1379B
5067-5383	Vacuum chamber replacement kit Includes two 5067-4798 and one 5067-5380	For G1379A, G1379B



1200 Series Vacuum Degasser



0100-1710



8710-1930



5063-6598



5063-6599

TECH TIP | Vacuum Degasser Care

- ➔ To generally clean the vacuum degasser tubing, flush the system with isopropanol
- ➔ Flush the degasser with water after using buffers
- ➔ Speed solvent changes by drawing solvent through the degasser and tubing with syringe adapter kit

Agilent Injection Valves

AGILENT INJECTION VALVES							
VALVE DESCRIPTION	USE WITH	VALVE	RHEBUILD KIT	ROTOR SEAL MATERIAL	ROTOR SEAL	STATOR	STATOR FACE
2 position/6 port injection valve, 400 bar	G1313A, G1329A, G1367A/B, 1120	0101-0921	0101-1257	Vespel	0100-1853	0100-1850	0100-1851
				Tefzel	0100-1849		
				PEEK	0100-2231		
2 position/6 port micro injection valve, 400 bar	G1377A, G1389A	0101-1050	—	Vespel	0100-2088	0100-2089	—
2 position/6 port injection valve, 600 bar	G1329B, G1367C SL, G1367D SL Plus, G1367E, 1220 Infinity LC	0101-1422	—	PEEK	0101-1416	0101-1417	—
2 position/6 port MBB injection valve, 400 bar	G2260A	—	0101-1268	Peek	0101-1268*	0100-2195	—
2 position/6 port ultra high pressure valve, 1200 bar	G4226A	—	—	Vespel	5068-0007	5068-0006	—

*Includes seal and stator face

0101-0921



0101-1417



5068-0006

