


Rheodyne RheFlex® Fingertight Fittings

The Rheodyne RheFlex Precision Two-Piece PEEK Fittings sets provide inert, biocompatible connections for instrumentation. These fittings have a reliable, time tested design. Also offered in this product line is the 6000-078 fitting, designed to connect 1/8" OD tubing into our manual preparative-scale injection valves.

RHEODYNE RHEFLEX FINGERTIGHT FITTINGS

PART NO	DESCRIPTION	TUBING OD	PORT	PRESSURE RATING	HEAD STYLE	MATERIAL	QTY
RHEFLEX TWO-PIECE FITTINGS (INCLUDES FERRULES)							
6000-078	Fitting	1/8"	5/16-24 Coned	5,000 PSI	5/16" Hex	PEEK, Nat	ea
6000-079	Repl Ferrules	1/8"	5/16-24 Coned	7,000 PSI	ChromTRAC knob	PEEK, Nat	5/pk
6000-054	Fitting	1/16"	10-32 Coned	7,000 PSI	ChromTRAC knob	PEEK, Nat	5/pk
6000-254	Fitting	1/16"	10-32 Coned	7,000 PSI	ChromTRAC knob	PEEK, Nat	10/pk
6000-255	Fitting, Short	1/16"	10-32 Coned	7,000 PSI	1/4" Hex	PEEK, Nat	10/pk
6000-051	Repl Ferrules	1/16"	10-32 Coned	7,000 PSI	ChromTRAC knob	PEEK, Nat	5/pk
6000-251	Repl Ferrules	1/16"	10-32 Coned	7,000 PSI	ChromTRAC knob	PEEK, Nat	10/pk
RHEFLEX ONE-PIECE FITTINGS							
6000-282	Fitting	1/16"	10-32 Coned	5,000 PSI	ChromTRAC knob	PEEK, Nat	10/pk
6000-360	Fitting	1/32"	M4 Coned	5,000 PSI	1/4" Hex	PEEK, Nat	10/pk



Upchurch Scientific® Fingertight Fittings

UPCHURCH SCIENTIFIC ONE-PIECE FINGERTIGHT FITTINGS

PART NO	DESCRIPTION	TUBING OD	PORT	PRESSURE RATING	INCLUDES FERRULE	QTY
F-100	PCTFE, Red, One-Piece Fitting	1/16"	10-32 Coned	4,000 PSI	N/A	ea
F-120	PEEK One-Piece Fitting	1/16"	10-32 Coned	5,000 PSI	N/A	ea
F-130	PEEK One-Piece Fitting, Long	1/16"	10-32 Coned	5,000 PSI	N/A	ea
M-645	PEEK/PCTFE One-Piece Fitting	1/32"	6-40 Coned	1,750 – 3250 PSI	N/A	ea
P-100	PCTFE One-Piece Fitting	1/8"	1/4-28 Coned	1,000 PSI	N/A	ea

UPCHURCH SCIENTIFIC TWO-PIECE FINGERTIGHT FITTINGS

F-200	Delrin Nut, Single Wing	1/16"	10-32 Coned	6,000 PSI	F-142 (PEEK)	ea
F-300	PEEK Nut, Double Wing	1/16"	10-32 Coned	6,000 PSI	F-142 (PEEK)	ea
F-330	PEEK Nut, Standard Knurl, Long	1/16"	10-32 Coned	6,000 PSI	F-142 (PEEK)	ea
F-330	PEEK Nut, Standard Knurl	1/16"	10-32 Coned	6,000 PSI	F-142 (PEEK)	ea

FERRULE OPTIONS FOR TWO-PIECE FINGERTIGHT FITTINGS

F-113	PEEK Natural Ferrule	1/32"	10-32 Coned	6,000 PSI	—	ea
F-142	PEEK Natural Ferrule	1/16"	10-32 Coned	6,000 PSI	—	ea
F-142N	ETFE Natural Ferrule	1/16"	10-32 Coned	4,000 PSI	—	ea
F-148	PCTFE Natural Ferrule	190 µm	10-32 Coned	6,000 PSI	—	ea
F-151	PCTFE Natural Ferrule	360 µm	10-32 Coned	6,000 PSI	—	ea
M-215	Conductive Ferrule	360 µm	10-32 Coned	1,500 PSI	—	ea

PLUGS

P-551	PEEK Natural Plug	—	10-32 Coned	—	N/A	ea
U-467R	Delrin Red Plug	—	10-32 Coned	—	N/A	ea

